



# AER Series

## AC20 - AutoCPAP



Respro AC20 Auto Cpap AER Series with HumidAir Heated Humidifier is a top-of-the-line Auto Cpap machine. This US FDA device features a design that makes it easy to use and maintain.

It features advanced features such as AutoSet Therapy, SmartStart, and AutoRamp, making it an excellent choice for those looking for reliable and effective Auto Cpap therapy.

The HumidAir Heated Humidifier is designed to provide optimal humidity levels for comfortable and effective Auto Cpap therapy.

.The machine is also compatible with a variety of masks, making it an excellent option for those looking for a complete Auto CPAP solution.

We are the only company in India to give 3 years warranty on AC20 Auto Cpap.



# RESPRO

# AER Series

AC20 - AutoCPAP

## Features

Modes	CPAP, Automatic CPAP (APAP)
CPAP Pressure range	4–20 cm H <sub>2</sub> O
Warranty	3 years
Includes:	



Nasal Mask



Heated Humidifier



Carry Bag



Hose Pipe



Free Gift

## Basic information

Model	AC20
Intended Use	OSA
User type	weight above 30kg

## Physical characteristics

Weight	1.35 kg
Size	238 x 178 x 128 mm
Screen	3.5 inches, color
Noise level (ISO80601-2-70)	≤ 28 dB

## Additional Parameters

Respiratory event detection	Leakage compensation
Ramp	Auto on/off
Auto Ramp	Mask fit
IPR (Intelligent pressure release)	Circuit compensation
Auto altitude adjustable	

## Humidifier

Integrated humidifier  
Auto humidity  
Preheating  
Heated tube

## Reports

Take out reports from the comfort of your Home.

Download Free software from our website on your laptop/ Desktop computer

SD card saves data upto 2-4 years  
Data storage in main unit  
Wireless transmission (WiFi, Bluetooth, GPRS)  
Multi-platform synchronous data sharing for patients and clinicians (Android, iOS, Web)

## Other Add-Ons

SD card	Free
Full face mask	Optional
Nasal pillow mask	Optional
Tube (15 mm)	Free
Tube (19 mm)	Optional
Heated tube (19 mm)	Optional
WiFi & Bluetooth module	Optional
Bluetooth SpO <sub>2</sub> watch	Optional
iMatrix intelligent cloud system	Optional



# RESPRO

Head Office: 1, National Mill Compound, Behind HP  
Petrol Pump, S.V.road, Dahisar East, Mumbai – 400064

☎ 9022120120 |

www.resproindia.com